

## Desalting of tryptic peptide mixtures using C18 Zip-tips

*This step is needed only for the direct analysis of tryptic peptide mixtures by MALDI mass spectrometry. If the sample will be prepared for MS analysis using a nano-HPLC or submitted to LC-ESI mass spectrometry, this step can be omitted.*

### Materials

Trifluoroacetic acid (Fluka)

Zip-tips, C18 (Millipore)

1. Dissolve dried peptide mixtures (obtained from steps 14 through 19) in 10  $\mu$ l of 5 % acetonitrile containing 0.1% trifluoroacetic acid.

For not completely dried peptide mixtures, add ultrapure water to make a volume of 8.5  $\mu$ l, and add 0.5  $\mu$ l acetonitrile and 1  $\mu$ l of 1% trifluoroacetic acid.

Check pH of resulting peptide mixture by using a pH paper; confirm and eventually adjust pH of samples to 2-4 using 1 % trifluoroacetic acid.

2. Prepare Zip-tip C18 (Millipore).  
Activate Zip-tip by rinsing 3 times with 10  $\mu$ l acetonitrile.  
Equilibrate Zip-tip 3 times with 10  $\mu$ l of 5 % acetonitrile containing 0.1 % trifluoroacetic acid.
3. Bind the peptides by aspirating and dispensing peptide solution; use up to 10 cycles.
4. Wash with 10  $\mu$ l of 5 % acetonitrile containing 0.1% trifluoroacetic acid.
5. Elute peptides with 10  $\mu$ l of 50 % acetonitrile containing 0.1 % trifluoroacetic acid. Carefully, aspirate and dispense eluant through Zip-tip without introducing air; repeat up to 10 times.

Alternatively, elute the digested peptide with 10  $\mu$ l of 20% acetonitrile containing 0.1% trifluoroacetic acid and with 10  $\mu$ l of 60% acetonitrile containing 0.1% trifluoroacetic acid. Combine these two eluants.

6. Dry eluant using a Speed-Vac centrifuge at 35 °C.
7. If necessary, store at -20 °C (< 24 hr) or -80 °C (for longer periods).
8. Prepare sample for MALDI mass spectrometric analysis:

Resolubilize the desalted peptide lyophilate with 1  $\mu$ l matrix solution (10  $\mu$ g/ $\mu$ l  $\alpha$ -cyano-4-hydroxy-cinnamic acid in 0.1% trifluoroacetic acid:acetonitrile, 1:1, v/v).

Load a half microliter of sample mixture on a MALDI target plate (with 144 spots for the AB 4700 instrument) and let dry at room temperature. The sample is now ready for MALDI mass spectrometry.

---

Oregon State University

Maier Lab – Biomolecular Mass Spec

Filename of document:DESALTING OF TRYPTIC PEPTIDE MIXTURES USING C18 ZIP TIPS